

Introduction

Welcome to The VisiVest™ System — a combination of the Bluetooth® enabled The Vest® System and The VisiView® Health Portal.

This connectivity guide will take you through the steps to begin securely transmitting your session data and logging into the VisiView® Health Portal so you and your health care team can work better together as a team.

Once connected, your therapy session details can be shared between you and your health care team, including:

- Duration – The time you've spent doing your therapy sessions.
- Settings – The frequency and pressure you're using.
- Cough Pause® Segments – Registers the number of program pauses taken during your therapy session.



Airway Clearance System

User Created Login

Account Number: _____

User Name: _____

Password: _____

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with health care providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health.** Hill-Rom's people, products and programs work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Bluetooth® is a registered trademark of Bluetooth SIG, Inc.

VisiView®, VisiVest™ and The Vest® are trademarks or registered trademarks of Hill-Rom Services PTE Ltd.

Hill-Rom, Enhancing Outcomes for Patients and Their Caregivers, and Cough Pause® are registered trademarks of Hill-Rom Services, Inc.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services PTE Ltd. ALL RIGHTS RESERVED.
192733 rev 4 16-JUN-2017 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

800-426-4224

www.respiratorycare.hill-rom.com

VisiView® Health Portal

Connectivity Guide for Patients and Caregivers



Enhancing outcomes for patients and their caregivers:

Hill-Rom

Enhancing outcomes for patients and their caregivers:

Hill-Rom

Transmitting data to the the VisiView® Health Portal

Note: VisiVest™ Air Pulse Generator needs to be plugged in/ powered on. The Data Hub should be within 30 feet of the air pulse generator. Select an unused outlet that is not on a switch to avoid accidentally turning it off. The Data Hub should be left plugged in.


- 1 Plug in your Data Hub to start the self-running activation process. The bar and home icons light up to indicate activity. This may take up to 5 minutes to complete.

SOLID YELLOWS



Plug in and Start Up
Automatic Activation

FLASHING YELLOWS



Initializing

SOLID GREEN/YELLOW




Ready for Pairing
with your Medical Device

- 2 The Data Hub pairs and collects medical data from your device.

FLASHING GREEN HOME



Note: For USB capable devices, first connect the Data Hub to your medical device using a USB cable



Data Transfer from Device in Progress (Begins Automatically)

- 3 The Data Hub sends collected data over a cellular network.

FLASHING GREEN BAR



Wireless Network Activity: Transmitting Data

FLASHING BLUE LIGHT



Successful Data Transmission

SOLID GREEN/YELLOW



Ready for Pairing with your Medical Device

Setting Up the User Portal

1. On your computer or mobile device open Firefox or Google Chrome internet browser. Please note: May not work with Internet Explorer. Go to the following website: <https://visiview.hill-rom.com>
2. Use the following user name and password for set up. (You can find this information on your training paperwork or on your packing slip.)

- **User Name:** Account # will be provided at time of training.
- **Password:** zip code, first 4 letters of your last name, date of birth (MMDDYYYY)

Example:

Zip Code: 55432

Last Name: Smith

Date of Birth: 10/15/1990

Password: 55432smit10151990



Note: If your name is less than 4 characters, enter all the characters of your last name.

3. You will be prompted to enter your email address, new password, and establish a security question/answer.
4. When you get logged in you will want to verify that your information matches what is on the screen. NOTE: If your information is inaccurate, please contact Hill-Rom customer service 1.800.426.4224.

How It Works



Patient uses The VisiVest™ System.



Information is available to view in the VisiView® Health Portal for you and your Health Care Team.

The Data Hub collects usage data from The VisiVest™ System and wirelessly delivers data. Information is stored on a secure server.